

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PRADC1/PAP21 Protein(C-6His)

产品货号: TP05920

第三版

别名:Protease-associated domain-containing protein 1;Protease-associated domain-containing protein of 21 kDa;Hpap21;C2orf7;PAP21;UNQ833/PRO1760

描述:PRADC1, also known as C2orf7 or PAP21, is short for Protease-associated domain-containing protein 1. It is a 188 aa. with a 21 aa. signal, and the 171 located Asn can be glycosylated. PRADC1 has two mutagenesis which are N121Q and N171Q. This protein is secreted and highly expressed in skeletal muscle, heart and liver. It is expressed at intermediate level in kidney.

配方:PBS

分子量:20.0 kDa

序列:His22-Trp188

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.