



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Vasorin/SLITL2/VASN Protein(C-6His)

产品货号: TP05924

第三版

别名:Vasorin; Protein Slit-Like 2; VASN; SLITL2

描述:Vasorin is a Type I membrane protein, which is predominantly expressed in vascular smooth muscle cells in a developmentally regulated pattern. The expression level of Vasorin can be down regulated during vessel repair after arterial injury. Vasorin binds to transforming growth factor beta (TGF- β) and attenuates TGF- β signaling in vitro.

配方:PBS

分子量:60.8 kDa

序列:Cys24-Pro575

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.