

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Angiopoietin-Like Protein 8/ANGPTL8/Betatrophin Protein(C-6His)

产品货号: TP05927

第三版

别名:Betatrophin;Angiopoietin-like protein 8;Lipasin;Angptl8

**描述:**The protein specifically promotes pancreatic beta cell proliferation and beta cell mass expansion, thereby improving glucose tolerance. It promotes pancreatic beta cell proliferation without insulin resistance. Also it acts as a blood lipid regulator by regulating serum triglyceride levels and possibly by promoting ANGPTL3 cleavage. It interacts with ANGPTL3. It predominantly expressed in liver and also expressed in adipose tissues. The ability of the protein to induce pancreatic beta cell proliferation is promising in diabetes therapy. Betatrophin treatment could supply or replace insulin injections by increasing the number of insulin-producing cells in diabetes.

**配方:**PBS

**分子量:**21.1 kDa

**序列:**Ala22-Ala198

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.