

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Angiotensinogen/Serpin A8/AGT Protein(C-6His)

产品货号: TP05938

第三版

别名:Angiotensinogen;Serpin A8; AGT and SERPINA8

描述:Angiotensinogen is involved in maintaining blood pressure and in the pathogenesis of essential hypertension and preeclampsia. Genetic variations in Angiotensinogen are associated with susceptibility to essential hypertension, and can cause renal tubular dysgenesis, a severe disorder of renal tubular development. Defects in the encoding gene have also been associated with non-familial structural atrial fibrillation, and inflammatory bowel disease.

配方:PBS

分子量:51.7 kDa

序列:Asp34-Ala485

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.