

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human dUTP Pyrophosphatase/dUTPase Protein

产品货号: TP05939

第三版

别名: Deoxyuridine 5'-Triphosphate Nucleotidohydrolase Mitochondrial; dUTPase; dUTP Pyrophosphatase; DUT

**描述:** Deoxyuridine 5'-Triphosphate Nucleotidohydrolase Mitochondrial (dUTPase) belongs to the dUTPase family. dUTPase exists as a homotrimer and is involved in nucleotide metabolism. dUTPase produces dUMP, the immediate precursor of thymidine nucleotides and it decreases the intracellular concentration of dUTP so that uracil cannot be incorporated into DNA. The dUTPase increase in PCR product yield, length and fidelity enables further down-stream applications. These effects make dUTPase useful in PCR fidelity and yield-sensitive applications. dUTPase is specific for dUTP and is critical for the fidelity of DNA replication and repair.

**配方:** PBS

**分子量:** 20.0 kDa

**序列:** Met 1-Asn164

**纯度:** > 95% by HPLC

**浓度:**

**内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.