

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Ankyrin Repeat and SOCS Box Protein 13/ASB13 Protein(N-6His)

产品货号: TP05944

第三版

别名: Ankyrin repeat and SOCS box protein 13; ASB13;

描述: Ankyrin repeat and SOCS box protein 13(ASB13) is a member of the ankyrin repeat and SOCS box-containing (ASB) family. ASB13 contain six ankyrin repeats sequence and a SOCS box domain. The SOCS box serves to couple suppressor of cytokine signalling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. ASB13 may be a substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

配方: PBS

分子量: 31.4 kDa

序列: Met 1-Asn278

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.