



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Plasminogen Activator Inhibitor 1/PAI-1/SERPIN E1 Protein(C-6His)

产品货号: TP05948

第三版

别名: Plasminogen activator inhibitor 1; Endothelial plasminogen activator inhibitor; Serpin E1; Mr1; Pai1; Planh1;

描述: Plasminogen activator inhibitor-1 (serpin E1) is a serine protease inhibitor which belongs to the serpin family. Serpin E1 acts as 'bait' for tissue plasminogen activator, urokinase, protein C and matriptase-3/TMPRSS7. Its rapid interaction with PLAT may function as a major control point in the regulation of fibrinolysis.

配方: PBS

分子量: 43.7 kDa

序列: Leu24-Ser402

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.