



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Plexin Domain-Containing Protein 2/PLXDC2 Protein(C-6His)

产品货号: TP05952

第三版

别名:Plexin domain-containing protein 2;Tumor endothelial marker 7-related protein;Plxdc2;Tem7r

描述:Mouse Plexin domain-containing protein 2(PLXDC2) is a single-pass type I membrane protein. The protein is expressed in tumor endothelium and in vessels of some normal tissues, such as the muscle and lung. PLXDC2 can interact with CTTN, and may play a role in tumor angiogenesis.

配方:PBS

分子量:49.0 kDa

序列:Glu31-Ala455

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.