

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Early Placenta Insulin-Like Peptide/INSL4/Placentin Protein(C-6His)

产品货号: TP05953

第三版

别名:Early Placenta Insulin-Like Peptide; EPIL; Insulin-Like Peptide 4; Placentin; INSL4

描述:Early Placenta Insulin-Like Peptide (INSL4) belongs to the insulin family. INSL4 is expressed in the early placental cytotrophoblast and syncytiotrophoblast. INSL4 is a secreted protein and a precursor that undergoes post-translational cleavage to produce 3 polypeptide chains, A-C, that form tertiary structures composed of either all three chains, or just the A and B chains. INSL4 plays an important role in the development of trophoblast and in the regulation of bone formation.

配方:PBS

分子量:13.6 kDa

序列:Ala26-Thr139

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.