

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Probable Serine Carboxypeptidase CPVL/VCP-Like Protein Protein(C-6His)

产品货号: TP05957

第三版

别名: Probable serine carboxypeptidase CPVL; Carboxypeptidase vitellogenic-like; Vitellogenic carboxypeptidase-like protein; VCP-like protein; CPVL

描述: Probable serine carboxypeptidase CPVL, also known as Carboxypeptidase, vitellogenic-like, Vitellogenic carboxypeptidase-like protein, is a member of the peptidase S10 family. It is expressed in macrophages but not in other leukocytes. And specifically, it is abundantly expressed in heart and kidney, also expressed in spleen, leukocytes, and placenta. This enzyme may be involved in the digestion of phagocytosed particles in the lysosome, and also participation in an inflammatory protease cascade, and trimming of peptides for antigen presentation.

配方: PBS

分子量: 52.9 kDa

序列: Leu23-Gly476

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.