



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Annexin A3/ANXA3 Protein

产品货号: TP05958

第三版

**别名:**35-alpha calcimedin;Annexin III;Annexin-3;Inositol 1;2-cyclic phosphate 2-phosphohydrolase;Lipocortin III;Placental anticoagulant protein III;PAP-III;ANX3

**描述:**Annexin A3 ( ANXA3 ) contains 4 annexin repeats and belongs to the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. It is abnormally expressed in fetuses of both IVF and ICSI, which may contribute to the increase risk of birth defects in these ART. ANXA3 is Inhibitor of phospholipase A2 and cleaves the cyclic bond of inositol 1,2-cyclic phosphate to form inositol 1-phosphate. it also also play a role in anti-coagulant properties.

**配方:**PBS

**分子量:**36.4 kDa

**序列:**Met 1-Asp323

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.