



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Progonadoliberein-2/GNRH2 Protein(C-6His)

产品货号: TP05967

第三版

别名:Progonadoliberein-2; Progonadoliberein II; Gonadoliberein-2; Gonadoliberein II; Gonadotropin-Releasing Hormone II; GnRH II; Luliberin II; Luteinizing Hormone-Releasing Hormone II; LH-RH II; GnRH-Associated Peptide 2; GnRH-Associated Peptide II; GNRH2

描述:Progonadoliberein-2, also known as Progonadoliberein II and GNRH2, belongs to the gonadotropin-releasing hormone family. GNRH2 is specially expressed in midbrain, at significantly higher levels outside the brain (up to 30-fold). GNRH2 can be cleaved into two chains, gonadoliberein-2 and GnRH-associated peptide 2. gonadoliberein-2 regulates reproduction in females by stimulating the secretion of both luteinizing- and follicle-stimulating hormones. The proproteins produce three transcript variants, but produce the same peptide hormone.

配方:PBS

分子量:11.5 kDa

序列:Gln24-Val120

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.