



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human V-Set and Transmembrane Domain-Containing 2A/VSTM2A Protein(C-6His)

产品货号: TP05968

第三版

别名:V-set and transmembrane domain-containing protein 2A;VSTM2

描述:VSTM2A, also called VSTM2, is short for V-set and transmembrane domain-containing protein 2A. It is a 236 aa. protein which can be departed into two parts. The first 24 aa. is a signal peptide, and the last is a V-set and transmembrane domain-containing protein 2A. There is a natural variant location which is E84K.

配方:PBS

分子量:25.5 kDa

序列:Ser24-Phe243

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.