



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human WW Domain-Binding Protein 1/WBP1 Protein(N-6His)

产品货号: TP05969

第三版

别名:WW Domain-Binding Protein 1; WBP-1; WBP1

描述:WW Domain-Binding Protein 1 (WBP1) is widely expressed in many tissues, but it is lowly expressed in the lung, placenta, kidney, and liver. WBP1 contains two WW-binding motifs: WW-binding 1 and WW-binding 2 that are involved in mediating protein-protein interactions through the binding of polyproline ligands. The WW-binding domain is composed of 38 to 40 semi-conserved amino acids shared by proteins with diverse functions including structural, regulatory, and signaling proteins. In addition, WBP1 also encodes a ligand of the WW domain of the Yes kinase-associated protein. This function has not been determined.

配方:PBS

分子量:12.6 kDa

序列:Gly170-Pro269

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.