



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human WW Domain-Binding Protein 2/WBP2 Protein(N-6His)

产品货号: TP05970

第三版

别名:WW Domain-Binding Protein 2; WBP-2; WBP2

描述:WW Domain-Binding Protein 2 (WBP2) is a ubiquitous protein that contains one GRAM domain. The WW domain is composed of 38 to 40 semi-conserved AA shared by proteins of diverse functions including structural, regulatory, and signaling proteins. The domain is participated in mediating protein-protein interactions. WBP2 binds to the WW domain of YAP1, WWP1 and WWP2. The WW-binding 1 motif of WBP2 mediates interaction with NEDD4. The function of this protein WBP2 has not been determined. Some researches demonstrate that WBP-2 also interacts with the thyroid-specific transcription factor Pax8.

配方:PBS

分子量:13.4 kDa

序列:Met 1-Ala100

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.