

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Apolipoprotein A-II/ApoA2 Protein(C-6His)

产品货号: TP05975

第三版

别名:Apolipoprotein A-II;Apo-AII;Apolipoprotein A2;Truncated apolipoprotein A-II;ProapoA-II;APOA2

描述:Apolipoprotein A-II(Apo-AII for short), also known as Apolipoprotein A2, is a secreted protein which belongs to the apolipoprotein A2 family. It exists as a disulfide-linked homodimer; and also can form a disulfide-linked heterodimer with APOD. APOA2 is the 2nd most abundant protein of the high density lipoprotein particles. This protein may stabilize HDL (high density lipoprotein) structure by its association with lipids, and affect the HDL metabolism. Defects in APOA2 gene might cause apolipoprotein A-II deficiency or hypercholesterolemia.

配方:PBS

分子量:9.7 kDa

序列:Gln24-Gln100

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.