

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Proprotein Convertase Subtilisin/Kexin Type 9/PCSK9 Protein(C-6His/HA/AVI)

产品货号: TP05982

第三版

别名:Proprotein Convertase Subtilisin/Kexin Type 9; Neural Apoptosis-Regulated Convertase 1; NARC-1; Proprotein Convertase 9; PC9; Subtilisin/Kexin-Like Protease PC9; PCSK9; NARC1

描述:Human Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) is a secretory subtilase belonging to the proteinase K subfamily. PCSK9 is synthesized as a soluble zymogen that undergoes autocatalytic intramolecular processing in the ER, the pro domain and mature chain secrete together through noncovalent interactions. PCSK9 binds with low-density lipoprotein receptor (LDLR) and plays a major regulatory role in cholesterol homeostasis. Inhibition of PCSK9 function by preventing PCSK9/LDLR interaction is currently being explored as a means of lowering cholesterol levels. PCSK9 also binds to apolipoprotein receptor 2 (ApoER2), and play a role in the neural development.

配方:PBS

分子量:14&59 kDa

序列:Gln31-Gln152&Ser153-Gln692

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.