

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Argininosuccinate Synthase/ASS1 Protein(N-6His)

产品货号: TP05999

第三版

别名: Argininosuccinate Synthase; Citrulline--Aspartate Ligase; ASS1; ASS

描述: Argininosuccinate Synthase (ASS1) is an urea cycle enzyme with a tetrameric structure composed of identical subunits. ASS1 is involved in the synthesis of arginine and catalyzes that condensation of citrulline and aspartate to argininosuccinate using ATP. ASS1 is important to the urea cycle as it catalyzes the important second last step in the arginine biosynthetic pathway. ASS1 mainly expressed in periportal hepatocytes, but also in most other body tissues. A deficiency of ASS1 causes citrullinemia (CTLN1), an autosomal recessive disease which is characterized by severe vomiting spells and mental retardation.

配方: PBS

分子量: 42.8 kDa

序列: Met 1-Lys412

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.