



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Protachykinin-1/TAC1 Protein(N-6His)

产品货号: TP06000

第三版

别名:Protachykinin-1;PPT;TAC1

描述:Protachykinin-1(TAC1) is a secreted protein and belongs to the tachykinin family. TAC1 is encoded by the TAC1 gene. This gene encodes four products of the tachykinin peptide hormone family, substance P and neurokinin A, as well as the related peptides, neuropeptide K and neuropeptide gamma. These hormones are thought to function as neurotransmitters which interact with nerve receptors and smooth muscle cells. They are known to induce behavioral responses and function as vasodilators and secretagogues.

配方:PBS

分子量:15.2 kDa

序列:Glu20-Arg129

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.