



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Zinc Finger Protein 70/ZNF70 Protein(N-6His)

产品货号: TP06004

第三版

别名: Zinc Finger Protein 70; Zinc Finger Protein N27C7-1; ZNF70

描述: Zinc Finger Protein 70 (ZNF70) is a member of the krueppel C2H2-type zinc-finger protein family. ZNF70 is localized to the cell nucleus and contains eleven C2H2-type zinc fingers. ZNF70 has sequence-specific DNA binding transcription factor activity, and it may be involved in transcriptional regulation. In addition, ZNF70 has the DNA-binding and metal-binding activity, which can bind zinc ion.

配方: PBS

分子量: 53.0 kDa

序列: Met 1-Leu446

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.