



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Protein CutA Protein(C-6His)

产品货号: TP06007

第三版

别名:Protein CutA; Acetylcholinesterase-Associated Protein; Brain Acetylcholinesterase Putative Membrane Anchor; CUTA; ACHAP; C6orf82

描述:Protein CutA (CUTA) possesses a signal peptide and is widely expressed in brain. CUTA may forms part of a complex of membrane proteins attached to acetylcholinesterase (AChE). CUTA takes part in cellular tolerance to a broad range of divalent cations other than copper. Alternate transcriptional splice variants, both protein-coding and non-protein-coding, have been found.

配方:PBS

分子量:17.9 kDa

序列:Met 1-Pro156

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.