

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Protein Disulfide-Isomerase A3/PDIA3 Protein(C-6His)

产品货号: TP06010

第三版

别名:PDIA3 protein;Protein disulfide-isomerase A3;58 kDa glucose-regulated protein;Disulfide isomerase ER-60;ER protein 57

描述:PDIA3 protein is also known as Protein disulfide-isomerase A3. It is a protein that in humans is encoded by the PDIA3 gene. PDIA3 is an enzyme that belongs to the endoplasmic reticulum and interacts with lectin chaperones calreticulin and calnexin to modulate folding of newly synthesized glycoproteins. PDIA3 interacts with thiazide-sensitive sodium-chloride cotransporter in the kidney and is induced by glucose deprivation. PDIA3 is part of the major histocompatibility complex (MHC) class I peptide-loading complex (TAP1), which is important for formation of the final antigen conformation and export from the endoplasmic reticulum to the cell surface.

配方:PBS

分子量:55.3 kDa

序列:Ser25-Leu505

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.