

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human ZW10 Interactor/ZWINT Protein(N-6His)

产品货号: TP06012

第三版

别名: ZW10 Interactor; ZW10-Interacting Protein 1; Zwint-1; ZWINT

描述: ZW10 Interactor is localized to the kinetochores from late Prophase to Anaphase and has a uniform distribution in the cytoplasm of Interphase cells. ZWINT interacts ZW10, MIS12 and NDC80 subunit of the NDC80 complex specifically during mitosis. ZWINT is a part of the MIS12 complex which is required for kinetochore formation and spindle checkpoint activity. In addition, ZWINT is required to target ZW10 to the kinetochore at prometaphase.

配方: PBS

分子量: 33.4 kDa

序列: Met 1-Pro277

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.