



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Protein Disulfide-Isomerase A5/PDIA5 Protein(C-6His)

产品货号: TP06016

第三版

别名:Protein disulfide-isomerase A5; Protein disulfide isomerase-related protein; PDIA5; PDIR

描述:Protein disulfide-isomerase A5 is a 519 amino acids protein that belongs to the protein disulfide isomerase family. It contains 3 thioredoxin domains. It can catalyze the rearrangement of -S-S- bonds in proteins.

配方:PBS

分子量:28.8 kDa

序列:Ser22-Leu262

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.