

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Syndecan-Binding Protein 1/Syntenin-1/SDCBP Protein(C-6His)

产品货号: TP06017

第三版

别名:Syntenin-1;Scaffold protein Pbp1;Syndecan-binding protein 1;Sdcbp

描述:SDCBP, also called syntenin-1, is short for Syndecan-binding protein 1. It is expressed by the gene Sdcbp. SDCBP seems to function as an adapter protein. In adherens junctions, the protein may function to couple syndecans to cytoskeletal proteins or signaling components. Meanwhile it seems to couple transcription factor SOX4 to the IL-5 receptor (IL5RA). SDCBP also play a role in vesicular trafficking, and is required for the targeting of TGFA to the cell surface in the early secretory pathway.

配方:PBS

分子量:33.4 kDa

序列:Ser2-Val299

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.