

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Leucine-Rich Repeat-Containing Protein 3B/LRRC3B Protein(C-6His)

产品货号: TP06036

第三版

别名:Leucine-Rich Repeat-Containing Protein 3B; Leucine-Rich Repeat Protein LRP15; LRRC3B; LRP15

描述:Leucine-Rich Repeat-Containing Protein 3B (LRRC3B) belongs to the LRRC3 family. LRRC3B is single-pass membrane protein and contains three leucine-rich repeats, one LRRCT domain, and one LRRNT domain. LRR-containing proteins, of which there are greater than 2,000, participate in many important processes, including plant and animal immunity, hormone-receptor interactions, cell adhesion, signal transduction, regulation of gene expression, and apoptosis. A number of microarray expression profiling studies on human cancers have shown that LRRC3B is down-regulated in gastric, breast, colon, testis, prostate and brain cancers, suggesting LRRC3B involvement in carcinogenesis

配方:PBS

分子量:20.0 kDa

序列:Cys34-Tyr204

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.