

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human BCAS2/DAM1 Protein(C-6His, N-T7 tag)

产品货号: TP06040

第三版

别名:Pre-mRNA-Splicing Factor SPF27; Breast Carcinoma-Amplified Sequence 2; DNA Amplified in Mammary Carcinoma 1 Protein; Spliceosome-Associated Protein SPF 27; BCAS2; DAM1

描述:Breast Carcinoma-Amplified Sequence 2 (BCAS2) is a member of the SPF27 family. BCAS2 is a nuclear protein and widely expressed in many tissues. BCAS2 is identified as being overexpressed in various breast cancer cell lines. BCAS2 is a component of the spliceosome, taking part in the removal of introns from mRNA precursors. BCAS2 interacts with estrogen receptor alpha and beta, thyroid hormone receptor beta, peroxisome proliferator-activated receptor gamma. BCAS2 functions as an ER co-activator and is capable of enhancing ER-mediated transcription.

配方:PBS

分子量:28.6 kDa

序列:Ala2-Phe225

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.