

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human ETS1/EWSR2/P54 Protein(N-6His)

产品货号: TP06044

第三版

别名:ETS1 Protein; Protein C-ets-1; V-ets Erythroblastosis Virus E26 Oncogene Homolog 1 (Avian); Isoform CRA_a; ETS1; hCG_1732038

描述:ETS1 Protein (ETS1) is a nuclear protein that belongs to the ETS family. Members of this family recognize the core consensus DNA sequence GGAA/T in target genes. Proteins function either as transcriptional activators or repressors of numerous genes. They are involved in stem cell development, cell senescence and death, and tumorigenesis. ETS1 is a transcription factor, containing one ETS DNA-binding domain and one PNT (pointed) domain. it has been shown to interact with TTRAP, UBE2I and Death Associated Protein 6.

配方:PBS

分子量:33.0 kDa

序列:Met 1-Thr 272

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.