



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Eukaryotic Translation Initiation Factor 4E/EIF4E Protein

产品货号: TP06052

第三版

别名:Eukaryotic translation initiation factor 4E;eIF-4E;eIF-4F 25 kDa subunit;mRNA cap-binding protein;EIF4E;EIF4EL1;EIF4F

描述:Eukaryotic translation initiation factor 4E is a 217 amino acids protein that belongs to the eukaryotic initiation factor 4E family. eIF4F is a multi-subunit complex, the composition of which varies with external and internal environmental conditions. It is composed of at least EIF4A, EIF4E and EIF4G1/EIF4G3. EIF4E is also known to interact with other partners. It can recognize and bind the 7-methylguanosine-containing mRNA cap during an early step in the initiation of protein synthesis and facilitates ribosome binding by inducing the unwinding of the mRNAs secondary structures.

配方:PBS

分子量:25.1 kDa

序列:Met 1-Val217

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.