

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Transglutaminase 2/TGM2 Protein(C-6His)

产品货号: TP06056

第三版

别名:Protein-glutamine gamma-glutamyltransferase 2; Tgm2; Tissue transglutaminase; Transglutaminase C; TGase-2

描述:Protein-glutamine gamma-glutamyltransferase 2 (TGM2) is a 78-kDa, calcium dependent enzyme,It belongs to the transglutaminase superfamily and transglutaminase family. The protein encoded by this TGM2 gene acts as a monomer, is induced by retinoic acid, and appears to be involved in apoptosis. TGM2 is the autoantigen implicated in celiac disease. Two transcript variants encoding different isoforms have been found for this gene.

配方:PBS

分子量:78.0 kDa

序列:Ala2-Ala686

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.