

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Retbindin/RTBDN Protein(C-6His)

产品货号: TP06080

第三版

别名:Retbindin; RTBDN

描述:Human Retbindin is a 229 amino acid secreted protein that belongs to the folate receptor family. The gene that encodes retbindin exists as two alternatively spliced isoforms. Retbindin is first expressed in retina. It may play a role in binding retinoids and other carotenoids as it shares homology with;riboflavin binding proteins. RTBDN gene was first identified in a study of human eye tissues.

配方:PBS

分子量:22.3 kDa

序列:Ser31-Pro229

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.