

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Fibroblast Growth Factor 2/FGF-2/FGFb (Met134-Ser288) Protein

产品货号: TP06094

第三版

别名:Fibroblast growth factor 2;FGF-2;Basic fibroblast growth factor;Bfgf;Heparin-binding growth factor 2;HBGF-2;FGF2;FGFB

描述:Fibroblast growth factor 2(FGF2) is a secreted protein and belongs to the heparin-binding growth factors family. FGF2 is produced by epithelial, tumor and other cell types. It involved in developmental processes and regulates differentiation, proliferation, and migration, FGF2 is a critical factor for growing embryonic stem cells in culture without inducing differentiation. FGF2 has a high affinity for heparan sulfate and binding is a step in the FGF basic activation of FGFR tyrosine kinase.

配方:PBS

分子量:17.2 kDa

序列:Met 134-Ser288

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.