



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Lysozyme G-Like Protein 2/LYG2 Protein(C-6His)

产品货号: TP06096

第三版

别名:Lysozyme G-Like Protein 2; LYG2; LYGH

描述:Lysozyme G-Like Protein 2 (LYG2) is a secreted protein that belongs to the glycosyl hydrolase 23 family. LYG2 contains one SLT domain, one protein domain present in bacterial lytic transglycosylase (SLT) and in eukaryotic lysozymes (GEWL). SLT domain catalyzes the cleavage of the β -1,4-glycosidic bond between N-acetylmuramic acid (MurNAc) and N-acetylglucosamine (GlcNAc). LYG2 has hydrolase activity which acting on glycosyl bonds, and possess lysozyme activity.

配方:PBS

分子量:22.6 kDa

序列:Ser20-Phe212

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.