



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human MANSC Domain-Containing Protein 1/MANSC1 Protein(C-6His)

产品货号: TP06119

第三版

别名:MANSC Domain-Containing Protein 1; Loss of Heterozygosity 12 Chromosomal Region 3 Protein; MANSC1; LOH12CR3

描述:MANSC domain-containing protein 1(MANSC1) is encoded by MANSC1 gene. MANSC1 is a Single-pass type I membrane protein which contains 1 MANSC domain. It is widely expressed in many tissues and mainly located in membrane.

配方:PBS

分子量:39.8 kDa

序列:Gln27-Leu385

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.