

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Fumarase/FH Protein

产品货号: TP06125

第三版

别名:Fumarate Hydratase Mitochondrial; Fumarase; FH

描述:Fumarase is an enzyme that catalyze the reversible hydration/dehydration of fumarate to S-malate and is involved in the tricarboxylic acid or Krebs cycle. Fumarase exists in both form, cytosolic form and N-terminal extend mitochondrial form. The N-terminal extended form is targeted to the mitochondrion, where the removal of the extension is the same form as in the cytoplasm. Fumarase is thought to act as a tumor suppressor, which deficiency can lead to progressive encephalopathy, cerebral atrophy and development delay.

配方:PBS

分子量:50.2 kDa

序列:Ala45-Lys510

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.