

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human MEK-Binding Protein 1/MP1/MAPKSP1 Protein(N-6His)

产品货号: TP06126

第三版

别名:Ragulator Complex Protein LAMTOR3; Late Endosomal/Lysosomal Adaptor and MAPK and MTOR Activator 3; MEK-Binding Partner 1; Mp1; Mitogen-Activated Protein Kinase Kinase 1-Interacting Protein 1; Mitogen-Activated Protein Kinase Scaffold Protein 1; LAMTOR3; MAP2K1IP1; MAPKSP1

描述:Mitogen-Activated Protein Kinase Scaffold Protein 1 (MAPKSP1) was identified as an interacting protein that belongs to the LAMTOR3 family. MAPKSP1 restricted to late endosomes by the mitogen-activated protein-binding protein-interacting protein, and binds specifically to MAP kinase kinase MAP2K1/MEK1, MAP kinase MAPK3/ERK1, and MAP kinase MAPK1/ERK2. MAPKSP1 interacts with MAP2K1/MEK1 and MAPK2 and enhances the activation of MAPK2, and thus is thought to function as an adaptor to enhance the efficiency of the MAP kinase cascade.

配方:PBS

分子量:15.8 kDa

序列:Met 1-Val124

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.