

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human GABA Receptor-Associated Protein-Like 1/GABARAPL1 Protein(N-6His)

产品货号: TP06130

第三版

别名:Gamma-Aminobutyric Acid Receptor-Associated Protein-Like 1; Early Estrogen-Regulated Protein; GABA(A) Receptor-Associated Protein-Like 1; Glandular Epithelial Cell Protein 1; GEC-1; GABARAPL1; GEC1

描述:Gamma-Aminobutyric Acid Receptor-Associated Protein-Like 1 (GABARAPL1) is a cytoplasmic protein that belongs to the MAP1 LC3 family. GABARAPL1 is expressed at very high levels in the brain, heart, peripheral blood leukocytes, liver, kidney, placenta, and skeletal muscle. It can interact with GABRG2, OPRK1 and α -Tubulin. GABARAPL1 increases cell-surface expression of kappa-type opioid receptor through facilitating anterograde intracellular trafficking of the receptor.

配方:PBS

分子量:16.2 kDa

序列:Met 1-Lys117

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.