



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Serpin A3/alpha-1-Antichymotrypsin Protein(C-6His)

产品货号: TP06136

第三版

别名:Alpha-1-Antichymotrypsin; ACT; Cell Growth-Inhibiting Gene 24/25 Protein; Serpin A3; SERPINA3; AACT

描述:Serpin A3 belongs to the Serpin superfamily of serine protease inhibitors. Serpin A3 has been shown to inhibit some serine proteases, such as neutrophil cathepsin G and mast cell chymase. Serpin A3 is synthesized initially in the liver and secreted in plasma. Serpin A3 has been found in the amyloid plaques from the hippocampus of Alzheimer disease brains. In addition to, Serpin A3 is associated with liver disease and Parkinson disease and chronic obstructive pulmonary disease.

配方:PBS

分子量:46.3 kDa

序列:His24-Ala423

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.