



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Galectin-14/LGALS14 Protein(C-6His)

产品货号: TP06146

第三版

别名:Placental Protein 13-Like; Charcot-Leyden Crystal Protein 2; CLC2; Galectin-14; Gal-14; LGALS14; PPL13

描述:Galectin-14 is a member of the Galectin family of carbohydrate binding proteins. The members of Galectin family contain one or two carbohydrate recognition domains, which can bind β -Galactoside. LGALS14 is expressed intracellularly in placenta and eosinophils, and is released by eosinophils following allergen stimulation. LGALS14 may be involved in the development of allergic inflammation. Two alternatively spliced transcript variants encoding distinct isoforms have been observed.

配方:PBS

分子量:17.1 kDa

序列:Met 1-Asp139

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.