

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Serpin B12/SERPINB12 Protein(C-6His)

产品货号: TP06147

第三版

别名:SERPINB12;Serpin B12 ; Serine Protease Inhibitor-clade B12

描述:Serpin B12 is a member of the serpin family. Serpins are the largest and most diverse family of serine protease inhibitors. Most serpins are secreted and attain physiologic concentrations in the blood and extracellular fluids. Serpin B12 is expressed in many tissues, including brain, bone marrow, lymph node, heart, lung, liver, pancreas, testis, ovary, and intestine. Serpins are involved in a number of fundamental biological processes such as blood coagulation, complement activation, fibrinolysis, angiogenesis, inflammation and tumor suppression and are expressed in a cell-specific manner. SerpinB12 inhibits trypsin and plasmin, but not thrombin, coagulation factor Xa, or urokinase-type plasminogen activator.

配方:PBS

分子量:49.5 kDa

序列:Met 1-Pro425

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.