

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Serpin B3/SERPINB3/SCCA1 Protein(N-6His)

产品货号: TP06148

第三版

别名:Serpin B3; Protein T4-A; Squamous cell carcinoma antigen 1; SCCA-1;serine (or cysteine) proteinase inhibitor, clade B (ovalbumin), member 3; serpin peptidase inhibitor, clade B (ovalbumin), member 3; Squamous cell carcinoma antigen 1; T4-A;SCCA1

描述:Serpin B3, also known as squamous cell carcinoma antigen-1 (SCCA-1), is a member of the serpin superfamily of serine protease inhibitors. Serpin B3 belongs to the subgroup ovalbumin-related serpins which are involved in the regulation of apoptosis, inflammation, angiogenesis and embryogenesis. It may act as a papain-like cysteine protease inhibitor to modulate the host immune response against tumor cells. It also functions as an inhibitor of UV-induced apoptosis via suppression of the activity of c-Jun NH(2)-terminal kinase (JNK1).

配方:PBS

分子量:46.0 kDa

序列:Met 1-Pro390

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.