

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Cathepsin B Protein

产品货号: TP06149

第三版

别名:Cathepsin B;Ctsb;Cathepsin B1

描述:Cathepsin B (CatB) is an enzymatic protein belonging to the peptidase (or protease) families. It is the first described member of the family of lysosomal cysteine proteases. It is known to process a number of proteins, including pro and active caspases, prorenin and secretory leucoprotease inhibitor (SLPI). It is believed to participate in intracellular degradation and turnover of proteins. Cathepsin B may play a role in activation and inactivation of caspases, activation of renin and inactivation of SLPI, the key steps in apoptosis, angiotensin production, and progression of emphysema, respectively. Cathepsin B may also have a role in invasion and metastasis.

配方:PBS

分子量:36.4 kDa

序列:Met 1-Phe 339

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.