

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Serpin B6/Placental Thrombin Inhibitor Protein(N-Trx, 6His)

产品货号: TP06151

第三版

别名: Serpin B6; Cytoplasmic Antiproteinase; CAP; Peptidase Inhibitor 6; PI-6; Placental Thrombin Inhibitor; SERPINB6; PI6; PTI

描述: Serpin B6 belongs to the serpin family. Serpin B6 localizes to the cytoplasm. Serpin B6 is expressed in many tissues, abundantly by mast cells in different tissues and mastocytoma lesions. Serpin B6 may be involved in the regulation of serine proteinases present in the brain or extravasated from the blood. In addition, Serpin B6 may play an important role in the inner ear in the protection against leakage of lysosomal content during stress and loss of this protection results in cell death and sensorineural hearing loss.

配方: PBS

分子量: 60.7 kDa

序列: Met 1-Pro376

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.