



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human MOB Kinase Activator 1A/MOB1A Protein(C-6His)

产品货号: TP06156

第三版

别名: MOB Kinase Activator 1A; Mob1 Alpha; Mob1A; Mob1 Homolog 1B; Mps One Binder Kinase Activator-like 1B; MOB1A; C2orf6; MOB4B; MOBK1B; MOBKL1B

描述: MOB Kinase Activator 1A (MOB1A) belongs to the MOB1/phocein family, which acts as an activator of LATS1/2 in the Hippo signaling pathway, plays a key role in organ size control and tumor suppression by restricting proliferation and promoting apoptosis. MOBKL1B stimulates the kinase activity of STK38 and STK38L. MOBKL1B binds to and regulate downstream targets such as the NDR-family protein kinases and LATS1 kinase.

配方: PBS

分子量: 26.0 kDa

序列: Ser2-Arg216

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.