



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human MOB-Like Protein Phocein/MOB4/MOBKL3 Protein(N-6His)

产品货号: TP06157

第三版

别名: MOB-Like Protein Phocein; 2C4D; Class II mMOB1; Mob1 Homolog 3; Mob3; Mps One Binder Kinase Activator-Like 3; Preimplantation Protein 3; MOB4; MOB3; MOBKL3; PHOCN; PREI3

描述: MOB-Like Protein Phocein is a member of the MOB1/Phocein Family. MOB-Like Protein Phocein is associated with membranes and the Golgi stacks. It is present in the cytosol, where it behaves as a protein complex. It has been shown that MOB-Like Protein Phocein interacts with DNM1, EPS15 and Nucleoside Diphosphate Kinase. MOB-Like Protein Phocein is the major partner of Striatin Family members and may play a important role in membrane trafficking, specifically in membrane budding reactions.

配方: PBS

分子量: 28.2 kDa

序列: Met 1-Ala225

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.