

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human C-C Motif Chemokine 18/CCL18/PARC Protein(N-6His)

产品货号: TP06159

第三版

别名:C-C Motif Chemokine 18; Alternative Macrophage Activation-Associated CC Chemokine 1; AMAC-1; CC Chemokine PARC; Dendritic Cell Chemokine 1; DC-CK1; Macrophage Inflammatory Protein 4; MIP-4; Pulmonary and Activation-Regulated Chemokine; Small-Inducible Cytokine A18;CCL18; AMAC1; DCCK1; MIP4; PARC; SCYA18

描述:C-C Motif Chemokine 18 (CCL18) is secreted protein that belongs to the intercrine beta (chemokine CC) family. CCL18 is expressed at high levels in the lung, lymph nodes, placenta, bone marrow, and dendritic cells. CCL18 is a chemotactic factor that attracts lymphocytes but not monocytes or granulocytes. CCL18 may be involved in B-cell migration into B-cell follicles in lymph nodes. CCL18 attracts naive T-lymphocytes toward dendritic cells and activated macrophages in lymph nodes. It has chemotactic activity for naive T-cells, CD4+ and CD8+ T-cells and thus may play a role in both humoral and cell-mediated immunity responses.

配方:PBS

分子量:10.1 kDa

序列:Ala21-Ala89

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.