

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Serpin Kazal-7/SPINK7/ECG2 Protein(C-6His)

产品货号: TP06163

第三版

别名:SPINK7;ECRG-2;Esophagus cancer-related gene 2 protein;Serine protease inhibitor Kazal-type 7;ECG2

描述:Serine protease inhibitor Kazal-type 7(SPINK7) is a secreted protein, that in humans is encoded by the SPINK7 gene. SPINK7 contains 1 Kazal-like domain. SPINK7 is probably serine protease inhibitor. Recombinant human SPINK7 is a single, non-glycosylated polypeptide chain containing 74 amino acids and fused to His-tag at c-terminus.

配方:PBS

分子量:8.2 kDa

序列:Ser20-Cys85

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.