



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Motilin Protein(C-6His)

产品货号: TP06165

第三版

别名:Promotilin; Motilin-associated peptide; MLN

描述:Promotilin is a secreted protein that belongs to the motilin family. Human Promotilin cDNA encodes a 115 amino acid (aa) precursor that contains a 25 aa signal peptide and a 90 aa with motilin and motilin-associated peptide. It plays an important role in the regulation of interdigestive gastrointestinal motility and indirectly causes rhythmic contraction of duodenal and colonic smooth muscle.

配方:PBS

分子量:11.4 kDa

序列:Phe26-Lys115

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.