

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Myelin Protein P0/MPZ Protein(C-6His)

产品货号: TP06176

第三版

别名: Myelin Protein P0; Myelin Peripheral Protein; MPP; Myelin Protein Zero; MPZ

描述: Myelin Protein P0 (MPZ) is a single-pass type I membrane glycoprotein which belongs to the myelin P0 protein family. MPZ contains one Ig-like V-type (immunoglobulin-like) domain, absent in the central nervous system. MPZ is a major component of the myelin sheath in peripheral nerves. It is postulated that MPZ is a structural element in the formation and stabilisation of peripheral nerve myelin, holding its characteristic coil structure together by the interaction of its positively-charged domain with acidic lipids in the cytoplasmic face of the opposed bilayer, and by interaction between hydrophobic globular of adjacent extracellular domains. Defects in MPZ associated with Charcot-Marie-Tooth disease and Dejerine-Sottas disease.

配方: PBS

分子量: 15.2 kDa

序列: Ile30-Arg153

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.