

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Resistin Protein (Fc Tag)

产品货号: TP06217

第三版

别名:ADSF;FIZZ3;RETN1;RSTN;XCP1

描述:Resistin is an adipocytokine, which has been studied for its role in insulin resistance and recently in inflammation. The RETN and CAP1 polymorphisms and gene expression may be potential biomarkers for breast cancer risk. Resistin (RETN), recently found to be relevant to inflammation and inflammatory disorders.

配方:PBS

分子量:38.2 kDa

序列:Ser17-Pro108

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.